

## Guidelines


This quick reference guide explains how to amend a Guideline in Vision 3. For more information on Guideline design - see [Guideline help](#).

### Finding the name of your Guideline

To amend a guideline, you must know its Mnemonic (short name). If you do not know what it is called, follow the steps below:


1. From the **Toolbar**, select **Local Guideline Index** .

*Select Local Guideline Index*




2. From the **Local Guideline Index** screen, select **Maintain** , to switch into design mode.

*From the Local Guideline Index Select  
Maintain*



3. Right click on the Guideline you want to amend, select **Edit** and make a note of the name.

Right click on the  
Guideline and  
select Edit



4. Now. Select **Cancel** to close the Guideline.

5. Select **Maintain**  to come out of design mode.

## Amending your Guideline

1. From Consultation Manager, select a test patient.
2. From the Toolbar choose **Guidelines - Select Guideline:**

*Choose Guidelines - Select  
Guideline*



3. Type the name of your guideline in the **Mnemonic** box and select **View**



4. The guideline displays as if you were going to review, or add data. To update the guideline you must be in design mode, select **Maintain**



5. The guideline displays in design mode with the floating toolbar:

*Design mode*

**Care Planning Consultation**

**Assessment of Needs**  
All Medical History Records in Descending order Where READ\_CODE = "389"

**Care Plan**  
Last 2 Medical History Records Where READ\_CODE = "8CM"

**Examination**  
Last 1 Blood pressure Records  
Last 1 Height Records  
Last 1 Weight Records  
Last 1 Pulse (CVS/BP) Records  
Last 1 Pulse oximetry Records

Last 1 Smoking Records  
Last 1 Alcohol Records


Last 1 Medical History Records Where READ\_CODE = "68e"

Last 1 Medical History Records Where READ\_CODE = "173H.00" OR READ\_CODE = "173I.00" OR READ\_CODE = "173J.00" OR READ\_CODE = "173K.00" OR READ\_CODE = "173L.00"

631 173H.003MRC Breathlessness Scale: grade 1 173I.003MRC Breathlessness Scale: grade 2 173J.003MRC Breathlessness Scale: grade 3 173K.003MRC Breathlessness Scale: grade 4 173L.003MRC Breathlessness Scale: grade 5 1

**Referrals**  
631 8Hk1.003Referral to falls service 1 8H440003Referral to cardiology multidisciplinary team 1 8H71.003Referral to dietician 1 8H76.003Refer to dietician 1 8H1F.003Referral to health trainer 1 8HkX.003Referral to exercise referral programme 1 8H62.003Referral to G.P. 1 ZL62.003Referral to clinical nurse specialist 1 8Hr.003Referral to health care assistant

6. For the purpose of this quick reference guide, we are only going to cover how to amend a patient data line and change the properties. For more information on the other options available click [here](#).
7. To change patient data, for example, to display the last 3 Blood pressure records (instead of the last 1 Blood pressure). Right click on the appropriate row and select **Edit**.
8. The patient data screen displays:



9. To display the last 3 Blood pressure records, replace last **1** Records with **3**  and select **Apply** .

**Note** - To further restrict the display of data in the guideline select **Filter**  and choose your criteria. For example you may want to include a date range.


10. Select **OK** to save the changes.

11. To check that the data displays correctly, select **Maintain**  to come out of design mode (which also saves the changes) and enables you view the data in the guideline.


## Properties

Sometimes you want to change how data displays in your Guideline. For example, you may want data to display directly under a heading without having to click on it to expand it. To do this:

1. Make sure you are in design mode in your guideline by selecting **Maintain** .
2. Right click on the row and choose **Properties**:




3. You can either select **Initially Expanded** so that items expand and display or change the **Level** (it depends on where and what row you are selecting in your guideline).



Setting a row to **Level 1** usually displays the data on the guideline (without having to select a heading to view it).

*When viewing the guideline the data  
below the heading: Examination  
displays immediately*



**4.** After making the line property changes, select **Close**.

**5.** Select **Maintain**  to save the changes and exit design mode.

## Useful Links

Some other useful links include:

- **Exporting guidelines**  
[http://help.visionhealth.co.uk/DLM660/Consultation\\_Manager/index.htm#24933](http://help.visionhealth.co.uk/DLM660/Consultation_Manager/index.htm#24933)
- **Importing guidelines**  
[http://help.visionhealth.co.uk/DLM660/Consultation\\_Manager/index.htm#24882](http://help.visionhealth.co.uk/DLM660/Consultation_Manager/index.htm#24882)