

# DLM 350

Appointment Enhancements

Dispensing Label Enhancements

New Cardiovascular and Stroke Risk Calculators

Test Requesting Changes

Northern Ireland – CCG Referrals

Wales – My Health Online Release 0

VOS – Online Appointments Changes

England - EPS Release 2 &SSO Changes

Consultation Manager Changes

Registration/Registration Links Changes

Mail Manager Changes



### Table of Editions and Contents

Date	Version	Contents	Output
20/05/2011	001		
28/06/2011	002	SSO Changes Added	Docx & PDF
21/07/2011	003	Recalls and WCCG	PDF
25/08/2011	004	General Edit	PDF
05/09/11	005	Added details VOS changes to What's New.	Docx & PDF
02/11/2011	006	Branch option for Smartcard users	Docx & PDF

Copyright © INPS Ltd 2011

---

---

# Contents

<b>DLM 350</b>	<b>1</b>
<b>Summary of Changes</b>	<b>1</b>
Consultation Manager	1
Registration & Registration Links	3
Mail Manager	3
Control Panel - Pharmacist Prescribing Code	4
Install .Net Framework	4
England - Spine Message Failure Notification	4
England - SSO	4
<b>Consultation Manager Changes</b>	<b>5</b>
Display NHS/CHI/H+C Number	5
Prevention Display Changes	6
Outstanding Recalls for General History and Procedure, Specimen or Sample	7
England - EPS Release 2 – Phase 3	7
Controlled Drugs	7
View Rejected and Cancel Failed	7
Mark Prescriptions as Cancelled	7
Cancellation Error Information	8
Bulk Signing	8
Edit Information for Patient/Dispenser	9
Reminder - Printing Tokens for Previous Electronic Prescriptions	10
Creating a WCCG Referral from the Floating Toolbar	11
SCI Pathology - Scotland	11
Structured Data Area (SDA) Changes	12
Numeric Value Added	12
Cardiac Enzymes	13
Death Administration	14
Faecal Occult Blood	14
Knee X-Ray	14
Lung Function	15
Scoring Test Result	15
Total Protein	15
Tri-Iodothyronine (T3)	15
Upper Leg X-Ray	15
Absence of Condition	16
Contraception	16
Drug Allergy and Intolerance	17
Family History	18
Maternity	19
<b>Registration Changes</b>	<b>20</b>
Previous Address and Previous GP	20
<b>Registration Links Changes</b>	<b>21</b>
Re-acceptance	21

<b>Control Panel</b>	<b>22</b>
Pharmacist Prescribing Code	22
<b>Options - Install .Net Framework</b>	<b>22</b>
<b>England - Spine Message Failure Notification</b>	<b>23</b>
To Enable Background Queue Monitoring	23
Unsent Message Warning	23
<b>England - SSO</b>	<b>24</b>
New Unlock Session Message	24
Surgery Selection at Login	24

---

# DLM 350

---

## Summary of Changes

The following pages list changes to Vision functionality in DLM 350. The following projects are also part of DLM 350 but are documented separately due to their size:

- Advanced Label Design
- New Stroke Risk and CVD Risk Calculators
- Electronic Test Requesting (now compatible with Technidata)
- Appointments Enhancements
- England - Full Electronic Prescribing Service Release 2 Functionality
- Wales – My Health Online
- VOS – Online Appointments (Changes to VO S Configuration and Disabling Account in Registration) See the updated documents at [http://www.inps4.co.uk/my\\_vision/vua/vos](http://www.inps4.co.uk/my_vision/vua/vos)
- Northern Ireland - Electronic Referrals

User guides for these projects can be found on the User Assistance DLM page at [http://www.inps4.co.uk/my\\_vision/vua/dlm/index.html](http://www.inps4.co.uk/my_vision/vua/dlm/index.html)

## Consultation Manager

### Display NHS/CHI/H+C Number

- There is a new option to display NHS, CHI or H+C number in the Title Format section of Consultation – Options – Setup – General. Once this is selected, the relevant patient number displays in the Consultation Manager title bar.

See “[Display NHS/CHI/H+C Number](#)” on page 5.

### Prevention Display

- **Target Blood Pressure** -In response to user requests the Target Blood Pressure is now visible alongside Blood Pressure on Summary - Prevention Display.
- **Peak Flow** and **Pulse Rate** are now displayed on the Prevention Display screen.

See “[Prevention Display Changes](#)” on page 6.

## Recall Change

- From DLM 350, when a recall is entered with Read code which is not mapped to any SDA (mapped to General History by default), it now auto-completes any previous recall with exactly the same code.

See "[Outstanding Recalls for General History and Procedure, Specimen or Sample](#)" on page 7.

## Label Trace

- The Label Trace icon has now been changed to .

## England - EPS Release 2 – Phase 3

- **Controlled Drug Changes** - All users can now see in the Options-Setup screen from the front screen of Vision, if schedule 2 & 3 and private prescribing are enabled for EPS 2.
  - **Cancellation changes**
    - You can now view rejected and cancel failed prescriptions in the bulk signing view - rejections tab.
    - You can now access the cancellation reason, dispenser and dispenser telephone number from Prescription Manager.
  - **Bulk Signing changes**
    - F9 to Sign** - You can now press F9 to sign selected prescription(s) in the workflow list.
    - Urgent Prescriptions** - Prescriptions marked as urgent are now displayed at the top of the selected patient's drug list in the bulk signing screen.
    - Prescriptions displayed oldest first** - The order of prescriptions in the patient drugs to sign in the bulk signing screen is now ordered oldest prescription first rather than alphabetical order.
    - Edit Additional Instructions** – You can now select Edit Additional Instructions from the bulk signing screen after a prescription has been sent to change information for the dispenser.
- See "[England - EPS Release 2 – Phase 3](#)" on page 7.

## Scotland – Scottish Cervical Call Recall Service (SCCRS)

- **SCCRS Recall Advice** - SCCRS Recall Advice - A "Refer for Colposcopy" code has been added to the current SCCRS recall advice codes.

## Scotland - SCI Pathology

- The address for the Highlands SCI Store has been updated to: <https://scistore1.nhsh.scot.nhs.uk/storews/store41/scistoreservices.aspx>

## Welsh Clinical Communications Gateway

- **WCCG referral from floating toolbar** – You can now trigger a referral by dragging and dropping your history entry onto the referral icon in the floating toolbar. See [“Creating a WCCG Referral from the Floating Toolbar”](#) on page 11

- 

## Structured Data Area (SDA) Changes

- We have also changed, added and removed some Read codes and data entry options from existing SDA pick lists.

See [“Structured Data Area \(SDA\) Changes”](#) on page 12.

## Registration & Registration Links

### Previous address and previous GP

- Registration - For both Immigrant and Ex-service Acceptance Types, you now must complete the Previous Address and Previous GP details.

See [“Previous Address and Previous GP”](#) on page 20.

- Registration Links - Previous address and previous GP are now mandatory for re-acceptance transactions.

See [“Re-acceptance”](#) on page 21.

### Scotland – Chronic Medication Service

- **Pharmacy Telephone Number** – In the Registration Consent tab, the pharmacy telephone number now displays alongside the pharmacy name for CMS registered patients.

## Mail Manager

### New Units of Measure

The following units of measure are now available in the assign valid units list:

- UgFEU/ml
- GPL U/ml
- GPM U/ml
- mmol/mol Hb
- nmol/gHb/h
- ml/min/1.73m
- ml/min/1.73m\*2

## Control Panel - Pharmacist Prescribing Code

As per changes to pharmacy regulation, Vision now accepts the new GPhC prescribing numbers which are seven numeric digits starting with a 2 (eg 2123456). You are encouraged to enter this in Control Panel – File Maintenance as prescriptions prescribed by pharmacists with the old number may be rejected by the PPA.

([www.pharmacyregulation.org/ourpublicfocus/changestopharmacyregulation/index.aspx](http://www.pharmacyregulation.org/ourpublicfocus/changestopharmacyregulation/index.aspx))

See “[Pharmacist Prescribing Code](#)” on page 22.

## Install .Net Framework

There is a new selection in the Options menu for .net installation and version checking. This is a per workstation option for LAN practices only (.net installation for practices on the VES is managed by INPS).

See “[Options - Install .Net Framework](#)” on page 22.

## England - Spine Message Failure Notification

If you have background queue monitoring enabled and you have a problem with your Spine messaging, a system tray alert displays a warning message advising you to check your GPC machine and deal with any unsent EPS messages.

See “[England - Spine Message Failure Notification](#)” on page 23.

## England - SSO

- **New Unlock Session** Message - When attempting to unlock your workstation a change has been made to the displayed message.
- **Surgery Selection at Login** – For practices with branch surgeries, you are now prompted which surgery you are working from when logging in with your Smartcard.

See “[England - SSO](#)” on page 24.




## Consultation Manager Changes

### Display NHS/CHI/H+C Number

There is a new option to display NHS, CHI or H+C number in the Title Format section of Consultation – Options – Setup – General. Once this is selected, the relevant patient number displays in the Consultation Manager title bar.

1. From **Consultation Manager**, select **Consultation – Option – Setup**.
2. Select the **General** tab.




*Consultation Manager Setup – General tab*

3. In the Consultation Manager Title Format section, select **Include NHS/CHI/H+C number**.
4. Click **OK** to save and close.



*Consultation Manager Title bar with NHS No displayed*

## Prevention Display Changes

We have made the following additions to the Prevention Display Screen (accessed from Consultation Manager from  or Summary – Prevention Display).

### Target Blood Pressure


Target Blood Pressure is now displayed alongside Blood Pressure on Summary - Prevention Display.

Target Blood Pressure displays from the latest record of either:

- 246K. Target systolic blood pressure
- 246L. Target diastolic blood pressure

### Peak Flow and Pulse Rate Display

The latest Peak Flow and Pulse Rate recorded are now displayed on the Prevention Display.



The screenshot shows a software window titled "Prevention Display" with a menu bar containing icons for home, help, and mouse. The main content area is titled "Prevention Details" and contains the following information:

- Cytology:**
- Inflammation:**
- Contraception:**
- BP:** 30 November 2005 120 / 70
- Target BP:** 09 June 2010 120 / 80
- Pulse:** 09 June 2010
- Smoking:** 30 November 2005 Never smoked 0 cigs/day
- Alcohol:** 30 November 2005 Current drinker 7 units/week
- Diet:** 30 November 2005 <None>
- Exercise:** 30 November 2005 <None>
- Weight:** kgs BMI= Height: metres
- Weight loss from baseline:** Weight loss: kgs
- Peak Flow:** 09 June 2010 4.9 L/min
- Urinalysis:**
  - Glucose:
  - Protein:

Below the main details are two tables for allergies and intolerances:

Date	Read Term	Drug	Reaction	Severity

Date	Read term	Reaction	Severity	Certainty	Notes

Three yellow boxes with black arrows point to the following fields:

- Target BP** (points to the Target BP row)
- Pulse** (points to the Pulse row)
- Peak Flow** (points to the Peak Flow row)

*Prevention Display with Target Blood Pressure, Peak Flow and Pulse Rate highlighted*

## Outstanding Recalls for General History and Procedure, Specimen or Sample

Previously when adding a recall, if an outstanding recall existed within a Structured Data Area (SDA), the existing recall status was changed to Complete and a new recall added. The exception to this was when adding recalls which belonged to the General History and Procedure, Specimens and Sample SDA's. Existing recalls for such entries were not updated automatically and had to be manually changed.

From DLM 350, when a recall is entered with Read code which is not mapped to any SDA (mapped to General History by default), it now auto-completes any previous recall with exactly the same code.

## England - EPS Release 2 – Phase 3

### Controlled Drugs



The following changes have been implemented for controlled drugs:

- **All** users can now see in the Options-Setup screen from the front screen of Vision, if schedule 2 & 3 and private prescribing are enabled for EPS 2. Ministerial approval is needed before this option can be used.

### View Rejected and Cancel Failed

You can now view **rejected** and **cancel failed** prescriptions from Bulk Signing - Rejections tab. This takes you to the rejected tab in Prescription Manager.

### Mark Prescriptions as Cancelled

When a prescription has been marked as cancelled, the All Prescription icon for the user it was assigned to displays as  which means that there are prescriptions that have been cancelled in your list. Previously, the icon display remains as such regardless of the date filter, from DLM 350 you can now mark prescriptions as cancelled which removes alert  for that patient. Simply right click on the item in the Prescription Manager window.

## **Cancellation Error Information**

For Cancellation errors, you can now access the cancellation reason, dispenser and dispenser telephone number:


1. In Consultation Manager, go to Therapy.
2. Display the Prescription Manager screen by pressing F9 or click on .
3. Find the prescription you would like to know about using the various filters.
4. Right click and select Properties. From here you can see the reason for the prescription cancellation along with the dispenser name and telephone number.

## **Bulk Signing**


### **F9 to Sign Prescriptions**

You can now press F9 to sign selected prescription(s) in the workflow list.

### **Urgent Prescriptions**


Prescriptions marked as urgent are displayed in the bulk signing screen at the top of the selected patient's prescription items list with the following symbol .


### **Prescriptions Displayed Oldest First**

The order of prescriptions in the patient's prescription items list in the bulk signing screen is ordered oldest prescription first, rather than in alphabetical order. Unless there are urgent prescriptions which are at the top of the list with the following symbol .

## Edit Information for Patient/Dispenser


From the bulk signing screen, you can now right click and select **Edit Additional Instructions**. This allows you to change information notes for the patient or dispenser:

1. In Consultation Manager, select the  **List All Prescriptions** icon (you do not have to have a patient selected to access this screen).
2. Click on the **Review Required** tab.
3. Select the required patient.
4. In the window that displays the prescription items right click.
5. Select **Edit Additional Instructions**.



*Bulk Sign - Edit Additional Instructions*

6. Change the notes as required and click OK.



*Bulk signing screen – patient selected – right click –Edit Additional Instructions*

### **Reminder - Printing Tokens for Previous Electronic Prescriptions**

To print a barcoded token for prescriptions which have been previously sent as electronic only:

1. In Consultation Manager, go to Therapy.
2. Display the Prescription Manager screen by pressing F9 or click on .
5. Find the prescription you would like to know about using the various filters. And tick the box to select it.
3. Now press Print.
4. A token is printed and no message is resent. The drug item is now displayed as electronic + token.

## Creating a WCCG Referral from the Floating Toolbar


You can trigger WCCG from the floating toolbar referral icon. This speeds up the referring process:

1. From **Consultation Manager**, with the patient selected and a consultation open, click and hold on the History entry you wish to refer.
2. The floating toolbar is displayed, point to the bowtie icon and release the mouse button, this opens the **Referral – Add** screen.



*Floating Toolbar - Make a referral for this item*

3. Complete the referral details and then click on the envelope icon to create a referral letter if required. Click on OK to save and close the referral.
4. You are prompted: "Do you want to add this referral to the WCCG?".



*Add WCCG message*

5. Click **Yes** to open WCCG, or **No** to return to Consultation Manager.

## SCI Pathology - Scotland

The address for the Highlands SCI Store has been updated to:  
<https://scistore1.nhsh.scot.nhs.uk/storews/store41/scistoreservices.asmx>

## Structured Data Area (SDA) Changes

### *Numeric Value Added*

The following test result Structured Data Areas (SDA's) now have a **Value** area for numeric results to be recorded, a **Unit of Measure** option and a **Normal Range** section. All three of the new areas are searchable.

SDA	Found within
B12 & Folate level	Biochemistry (Other)– Iron, B12 and Folate
Drug levels	Biochemistry - Drug Levels
Blood lipids	Biochemistry - Blood Lipids
Bone studies	Biochemistry - Calcium & Bone Studies
Liver Function tests	Biochemistry - Liver Function Tests
Liver enzymes	Biochemistry - Liver Function Tests
Female sex hormone profile	Biochemistry - Sex Hormones
Male sex hormone profile	Biochemistry - Sex Hormones
Gonadotrophin	Biochemistry - Sex Hormones
Thyroid function	Biochemistry - Thyroid Function Tests
Urea and Electrolytes	Biochemistry - Urea and Electrolytes
Serum electrolytes	Biochemistry - Urea and Electrolytes
Urine dipstick for glucose	Biochemistry - Urine Chemistry Tests
Urine dipstick for protein	Biochemistry - Urine Chemistry Tests
Urine dipstick for ketones	Biochemistry - Urine Chemistry Tests
Urine dipstick for blood	Biochemistry - Urine Chemistry Tests
Urine Dipstick for Nitrites	Biochemistry - Urine Chemistry Tests
Urine Leucocytes	Biochemistry - Urine Chemistry Tests
Toxicology	Biochemistry (Other) – Other Biochemistry Tests
Urinalysis - Glucose	Biochemistry (Routine)
Urinalysis - Protein	Biochemistry (Routine)
Spirometry	Diagnostic Tests - Lung Function
Blood group antibodies	Haematology - Blood Film & Group
Full blood count	Haematology - Full Blood Count
Clotting tests	Haematology - Other Haematology Results
RBC red blood cell shape	Haematology - Red Blood Cell Indices
RBC colour/staining	Haematology - Red Blood Cell Indices
RBC clumping	Haematology - Red Blood Cell Indices
RBC enzymes	Haematology - Red Blood Cell Indices
Ante natal blood tests	Maternity - Investigations
Urine test	Microbiology - Other Microbiology
Anti Nuclear antibodies	Serology & Immunology - Autoantibodies



SDA	Found within
Anti Mitochondrial antibodies	Serology & Immunology - Autoantibodies
Anti smooth muscle autoantibodies	Serology & Immunology - Autoantibodies
Parietal cell autoantibodies	Serology & Immunology - Autoantibodies
Adrenal autoantibodies	Serology & Immunology - Autoantibodies
DNA binding autoantibodies	Serology & Immunology - Autoantibodies
Other autoantibodies	Serology & Immunology - Autoantibodies
HLA Tissue Typing	Serology & Immunology - Autoantibodies
Syphilis test	Serology & Immunology – Infectious Disease Titres
SH-antigen (hepatitis B) test	Serology & Immunology – Infectious Disease Titres
HIV Test	Serology & Immunology – Infectious Disease Titres
Legionella antibody test	Serology & Immunology – Infectious Disease Titres
Infectious titres and antibodies	Serology & Immunology – Infectious Disease Titres
Hepatitis A Test	Serology & Immunology – Infectious Disease Titres
Hepatitis B Antibody	Serology & Immunology – Infectious Disease Titres
Viral studies	Serology & Immunology – Miscellaneous Serology
CAP (RAST)	Serology & Immunology – Miscellaneous Serology
Immunology Screening Tests	Serology & Immunology – Miscellaneous Serology
Serology	Serology & Immunology - Serology
Plasma electrophoresis	Serology & Immunology - Serology
LE cells	Serology & Immunology - Serology
Immuno-electrophoresis	Serology & Immunology - Serology

*Example of SDA with numeric result areas highlighted*

## Cardiac Enzymes

The following Read codes have been added to the available list of codes within the Cardiac Enzymes SDA:

- 44MC. Serum troponin T level
- 44ME. Plasma troponin I level
- 44MG. Serum troponin I level
- 44MH. Plasma troponin T level

### ***Death Administration***

The following inappropriate Read codes have been removed from the available list of codes within the Death Administration SDA:

- 94Z.. Death Administration NOS
- 94Z0. Preferred place of death
- 94Z1. Preferred place of death: home
- 94Z2. Preferred place of death: hospice
- 94Z3. Preferred place of death: community hospital
- 94Z4. Preferred place of death: hospital
- 94Z5. Preferred place of death: nursing home
- 94Z6. Preferred place of death: pt unable to express preference
- 94Z7. Preferred place of death: discussion not appropriate
- 94Z8. Preferred place of death: patient undecided
- 94Z9. Preferred place of death: learning disability unit
- 94ZA. Preferred place of death: mental health unit
- 94ZB. Preferred place of death: discussed with family

### ***Faecal Occult Blood***

The following Read codes have been added to the available list of codes within the Faecal Occult Blood SDA:

- 6867. BCSP faecal occult blood testing kit spoilt
- 6868. BCSP faecal occult blood test technical failure
- 6869. BCSP faecal occult blood test result unclear
- 686C. BCSP faecal occult blood test incomplete participation

### ***Knee X-Ray***

The following Read code has been added to the available list of codes within the Knee X-Ray SDA:

- 52A..13 Knee X-ray
- 52A8 Plain X-ray patella
- 52A80 Plain X-ray patella normal
- 52A81 Plain X-ray patella abnormal
- 52AB Stress X-ray knee

The following inappropriate Read code has been removed from the available list of codes within the Knee X-Ray SDA:

- 52A..14 Leg X-ray

## ***Lung Function***

In both **339a. Lung function before bronchodilation** and **339b. Lung function after bronchodilation** the value can now be recorded to two decimal places.

## ***Scoring Test Result***

The following Read code has been added to the available list of codes within the Scoring Test Result SDA:

- 687C. Malnutrition universal screening tool score

This makes the numeric result searchable for the first time.

## ***Total Protein***

The following inappropriate Read codes have been removed from the available list of codes within the Total Protein SDA:

- 44MC. Serum troponin T level
- 44ME. Plasma troponin I level
- 44MG. Serum troponin I level
- 44MH. Plasma troponin T level

## ***Tri-Iodothyronine (T3)***

The following Read code has been added to the available list of codes within the Tri-Iodothyronine (T3) SDA:

- 442f. Serum total T3 level

## ***Upper Leg X-Ray***


The following inappropriate Read codes have been removed from the available list of codes within the Upper Leg X-Ray SDA:

- 52A3. Plain X-ray hip joint
- 52A30 Plain X-ray hip joint normal
- 52A31 Plain X-ray hip joint abnormal
- 52A7. Plain X-ray knee
- 52A70 Plain X-ray knee normal
- 52A71 Plain X-ray knee abnormal
- 52A8. Plain X-ray patella
- 52A80 Plain X-ray patella normal
- 52A81 Plain X-ray patella abnormal
- 52A9. X-ray shaft of tibia/fibula
- 52A90 X-ray shaft of tibia/fibula normal
- 52A91 X-ray shaft of tibia/fibula abnormal
- 52A92 X-ray of tibia normal

- 52A93 X-ray of tibia abnormal
- 52A94 X-ray of fibula normal
- 52A95 X-ray of fibula abnormal
- 52AA. Plain X-ray ankle joint
- 52AA0 Plain X-ray ankle joint normal
- 52AA1 Plain X-ray ankle joint abnormal
- 52AB. Stress X-ray knee
- 52AC. Stress X-ray ankle joint

### **Absence of Condition**

The way Absence of Condition data appears on Guidelines and Reports has been changed to make the data clearer.



*Example of an Absence of Condition Guideline, data highlighted*

### **Contraception**

The following Read codes have been added to the available list of codes within the Contraception SDA:

- 7E09. Intrauterine contraceptive device procedure
- 7E090 Introduction of intrauterine contraceptive device
- 7E091 Replacement of intrauterine contraceptive device
- 7E092 Removal of intrauterine contraceptive device NEC
- 7E093 Removal of displaced intrauterine contraceptive device
- 7E094 Introduction of Mirena coil
- 7E095 Removal of Mirena coil
- 7E09y Other specified intrauterine contraceptive device
- 7E09z Intrauterine contraceptive device procedure NOS

## Drug Allergy and Intolerance

When adding a **Drug Allergy and Intolerance**, the **Severity** and **Certainty** no longer default. These fields must now be manually selected. This prevents inaccurate data being recorded in error.

The screenshot shows a form titled "Drug Allergy and Intolerance - Add". It includes fields for "Date of Recording" (23 February 2011), "Clinician" (Dr Jane Mars), and a checked "In Practice" box. There are dropdown menus for "Read Term for Allergy" (14L..00 H/O: drug allergy), "Drug", "Read Term for Reaction", and "Reaction Type" (Allergy). At the bottom, the "Severity" and "Certainty" dropdown menus are highlighted with a red rectangular box.

*Drug Allergies and Intolerance – Add with Severity and Certainty highlighted*

If you click **OK** without selecting a value for **Severity** and/or **Certainty**, the following messages display.

The dialog box is titled "Consultation Manager" and contains an information icon (i) on the left. The text reads: "Severity of allergy. An item must be selected from the list." On the right side, there are buttons for "Help" and "Notes >>". At the bottom center is an "OK" button.

*Severity error message*

The dialog box is titled "Consultation Manager" and contains an information icon (i) on the left. The text reads: "Certainty of allergy. An item must be selected from the list." On the right side, there are buttons for "Help" and "Notes >>". At the bottom center is an "OK" button.

*Certainty error message*

The following Read code has been added to the available list of codes within the Drug Allergy and Intolerance SDA:

- 8I7.. Treatment not tolerated

The following inappropriate Read codes have been removed from the available list of codes within the Allergy and Intolerance SDA:

- 8I2M. Spirometry reversibility testing contraindicated
- 8I2N. Drug dependence home detoxification contraindicated

### ***Family History***

The following Read codes can now be added through the Family History SDA:

- 122C0 No family history of asthma
- 124J. Family history of malignant melanoma
- 1254. Family history of Addison disease
- 1287.11 Family history of bipolar disorder
- 128A. FH: Psychotic illness
- 129C. Family history of subarachnoid haemorrhage
- 12CT. Family history of myocarditis
- 12CV. FH: Cardiovascular disease 1<sup>st</sup> degree male relative<55yrs
- 12CW. FH: Cardiovascular disease 1<sup>st</sup> degree female relative<65yrs
- 12D5. Family history of pulmonary emphysema
- 12H5. Family history of psoriasis with arthropathy
- 12W2. Family history of learning disability
- 12X.. Family history of substance misuse
- 12X0 Family history of alcohol misuse
- 12X1. Both parents misuse drugs
- 12X2. Paternal drug misuse
- 12Y.. Family history of sarcoidosis
- 12a.. Family history of alpha-1-antitrypsin deficiency

## Maternity

To save confusion between recording maternal details and infant details, three of the Maternity SDA's have been renamed:

- Delivery Details – Mother

**Delivery Details - Mother - Add**

Recall OK Cancel Help

Date of delivery: 16 March 2011 Date of recording: 16 March 2011 Private  In Practice  Clinician: Dr Clive Clover

Read Ierm: 63...00 Birth details Notes:

<Previous Next> Gestation - Maternity Outcome Skip to

*Delivery Details – Mother – Add screen*

- Pregnancy Outcome – Mother

**Pregnancy Outcome - Mother - Add**

Recall OK Cancel Help

Date: 16 March 2011 Clinician: Dr Clive Clover Private  In Practice  Read Ierm: 63...00 Birth details

Hospital: Discharged date: In Charge: Status:

Notes:

<Previous Next> Placenta Skip to

*Pregnancy Outcome – Mother – Add screen*

- Infant Details – Mother

**Infant details - Mother - Add**

Recall OK Cancel Help

Date: 16 March 2011 Clinician: Dr Clive Clover Private  In Practice  Read Ierm: 63...00 Birth details Date of birth:

Sex: <Not specified> Birth weight: Unit of measure: Kgs Apgar at 1 minute: Apgar at 5 minutes: Outcome: Live

Notes:

<Previous Next> Skip to


*Infant Details – Mother – Add screen*

---

## Registration Changes


### *Previous Address and Previous GP*

For both Immigrant and Ex-service Acceptance Types, you now must complete the Previous Address and Previous GP details. To prevent a registration being saved without this information, an \* is displayed to indicate mandatory information.



*Registration – General Previous Address\*and Previous GP\* highlighted*

If you do not enter a Previous Address and Previous GP, the following message is displayed "Previous address is mandatory for a Transfer-In". Please click OK and ensure a previous address and previous GP are entered.



*Registration – General error message*

Once the Previous Address and Previous GP are recorded the text reverts to blue.

---

**Tip** – 'Unknown' can be entered where the previous address or previous GP information is not available. However, your registration authority may reject this.

---



---

## Registration Links Changes

### *Re-acceptance*

The previous address and previous GP are now mandatory requirements for Re-acceptance transactions within Vision Registration Links.

When a re-acceptance is received from your Trading Partner (TP) following a rejection (this usually indicates the wrong TP was initially selected), and the previous address and/or the previous GP name have not been recorded:


1. From **Registration Links - Incoming Transactions**, select **Process**.
2. Click on the transaction that requires re-acceptance and select **Process**.
3. **Re-acceptance following the receipt of a Rejection** is displayed.
4. Change the **Responsible TP** as required from the available list.
5. Complete the **Previous Address**, if there is a previous entry, this option will be greyed out.
6. Complete the **Previous GP name**, if there is a previous entry, this option will be greyed out.

---

**Tip** – 'Unknown' can be entered where the previous address or previous GP information is not available

---

7. Enter GP Notes if required.
8. Click **OK** to save and submit.



*Example of Re-acceptance following the receipt of a Rejection screen with Previous Address greyed out*

If the Previous Address and Previous GP Name have been entered during initial registration there is no change.

---

## Control Panel

### Pharmacist Prescribing Code

As per changes to pharmacy regulation, Vision now accepts the new GPhC prescribing numbers which are seven numeric digits starting with a 2 (eg 2123456). You are encouraged to enter this in Control Panel – File Maintenance as prescriptions prescribed by pharmacists with the old number may be rejected by the PPA.

([www.pharmacyregulation.org/ourpublicfocus/changestopharmacyregulation/index.aspx](http://www.pharmacyregulation.org/ourpublicfocus/changestopharmacyregulation/index.aspx))

1. Go to Control Panel – File Maintenance – Staff.
2. Right click on the pharmacist and click Edit.
3. Go the Professional tab and replace the Prescribing code with the new GPhC prescribing number.
4. Press OK to save.

---

## Options - Install .Net Framework

For LAN practices, DLM 313 brought .net installation which is needed to run certain Vision functionality ie CVD Risk Calculators. From DLM 350, there is a new selection in the Options menu to Install .Net framework which allows you to view the versions of .net currently installed on the workstation. You can also select Install New Version to check for upgrades.

---

## England - Spine Message Failure Notification

From DLM 350, with Background Queue Monitoring Enabled, you are now warned if your messages are not sent from your Vision system to the Spine (ie if your N3 internet connection is down or there is a problem with your GPC machine).

### To Enable Background Queue Monitoring


1. From the main Vision screen, select **Options – Status**.
2. Click on the **Options** button.
3. Tick **Enable background message queue checks**.
4. Press **OK**.

### Unsent Message Warning

If you have a problem with your Spine messaging, a system tray alert displays the following:

*"There may be a problem with the messaging program. Spine services (PDS, EPS, CAB, etc) are unavailable. You are advised to check the GPC machine. Please also check any EPS prescriptions and take alternative actions as they may have not reached the spine."*

Please check that your GPC machine is working properly. When the GPC is working properly, you can then access Mail Manager and reprocess the messages to resend them.



---

**Note** - It is important you attend to the error message as soon as possible as your patient may be waiting for their EPS prescription at the pharmacy.

---

---


## England - SSO

### New Unlock Session Message

When your workstation is locked and you move your mouse or type on your keyboard, the Unlock Session screen displays, you are prompted:


*"Please remove and re-insert your smartcard to unlock the workstation. To end the current session instead, remove the smartcard and select 'End Session'.*

*Lock Reason: Vision maximum inactivity time reached."*



*Vision - Unlock Session*

To unlock the session you must remove and re-enter your Smartcard and passcode, you are then returned to Vision as it was when it was locked. Alternatively, you can click the **End Session** button to close Vision. A Vision warning is displayed concerning the loss of any unsaved data, click **Yes** to continue.



*Vision Warning*

### Surgery Selection at Login

For practices with branch surgeries, you are now prompted which surgery you are working from when logging in with your Smartcard.